West of Scotland Paediatric Gastroenterology, Hepatology and Nutrition Network



Training package for parents/carers

Replacement of Balloon Gastrostomy

Parent/Carer name	
Child's name	

	Date	Trainer signature	Parent/carer signature
Procedure discussed with parent/carer			
Procedure observed by parent/carer			
Procedure performed by parent/carer with supervision			
Procedure performed by parent/carer under observation by trainer			
Parent/carer is confident to practice procedure			
Trainer is confident that parent/carer is able to practise unsupervised			

Steps of procedure	Procedure discussed	Procedure observed	Procedure performed with supervision	Procedure performed under observation by trainer	Parent/carer is confident to practice procedure
How to acquire equipment and supplies					
Effective hand hygiene					
Collect equipment required and place on a clean dry surface					
Check replacement device is correct make, size, length and expiry date					
Check integrity of the balloon on the new gastrostomy to be inserted (using air)					
Explain procedure to child/young person					
Position the child/young person for changing the gastrostomy					
Deflate balloon of gastrostomy in situ					

Steps of procedure	Procedure discussed	Procedure observed	Procedure performed with supervision	Procedure performed under observation by trainer	Parent/carer is confident to practice procedure
Remove gastrostomy from the child/young persons stomach					
Clean the stoma area					
Apply lubricant gel to balloon end of gastrostomy and insert into stoma					
Inflate balloon with 3-5 mls of cooled, boiled water to secure in situ as instructed by your health care professional					
Attach extension set and aspirate with syringe to obtain 1-2mls of fluid from stomach					
Check pH reading of stomach fluid is 5.5 or less to ensure position of new gastrostomy in stomach					
Action to be taken if appropriate pH not obtained					

Steps of procedure	Procedure discussed	Procedure observed	Procedure performed with supervision	Procedure performed under observation by trainer	Parent/carer is confident to practice procedure
Disposal of waste					
Discuss possible problems including tube blockage					
How to seek help or advice					

Adapted by West of Scotland Paediatric Gastroenterology, Hepatology and Nutrition Managed Clinical Network (WoSPGHaN) with permission from NHS Lanarkshire Community Children's Nursing Team. If you would like more information, visit the website www.wospghan.scot.nhs.uk

WoSPGHaN

West of Scotland Gastroenterology, Hepatology and Nutrition Network www.wospghan.scot.nhs.uk

Published December 2012 Review December 2015 Version 1